



Mouse Anti-Human IL-34 monoclonal antibody, clone C-S45 (CABT-ZB387)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human IL-34
Target	IL34
Immunogen	Recombinant human IL-34
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-S45
Conjugate	Unconjugated
Applications	ELISA(Det) Recommended Working Concentration: ELISA: 0.05-0.5 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB385 - CABT-ZB387
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS
Preservative	None

Storage	Stable at +2-8°C for 12 months
----------------	--------------------------------

Ship	Wet ice
-------------	---------

GENE INFORMATION

Synonyms	IL34; interleukin 34; IL-34; C16orf77; interleukin-34; anti-IL-34
-----------------	---

Entrez Gene ID	146433
-----------------------	------------------------

UniProt ID	Q6ZMJ4
-------------------	------------------------