



# Mouse Anti-Human IL-6 monoclonal antibody, clone C-F9 [PE] (CABT-ZB396)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Recognizes both natural and recombinant human IL-6
<b>Target</b>	IL6
<b>Immunogen</b>	Recombinant human IL-6
<b>Isotype</b>	IgG1, κ
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	C-F9
<b>Conjugate</b>	PE
<b>Applications</b>	<p>ELISA(Cap), ELISPOT, FA, FC</p> <p>Recommended Working Concentration:</p> <p>ELISA: 0.5-2.5 µg/ml</p> <p>ELISpot: 5-10 µg/ml</p> <p>We recommend the following for sandwich ELISA (Capture - Detection):</p> <p>CABT-ZB394 - CABT-ZB397</p> <p>CABT-ZB396 - CABT-ZB397</p>
<b>Reconstitution</b>	Reconstitute with 1mL of deionized water.
<b>Format</b>	Purified, Lyophilised
<b>Size</b>	100 tests (1 mL)
<b>Buffer</b>	PBS

<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Stable at +2-8°C for six months after reconstitution. DO NOT FREEZE.
<b>Ship</b>	Wet ice

## GENE INFORMATION

<b>Synonyms</b>	IL6; interleukin 6; HGF; HSF; BSF2; IL-6; IFNB2; interleukin-6; CDF; BSF-2; IFN-beta-2; interferon beta-2; interleukin BSF-2; interferon; beta 2; hybridoma growth factor; CTL differentiation factor; B-cell stimulatory factor 2; B-cell differentiation factor; PF-04236921
<b>Entrez Gene ID</b>	<a href="#">3569</a>
<b>UniProt ID</b>	<a href="#">P05231</a>