



Mouse Anti-Human IL-6 monoclonal antibody, clone C-F5 (CABT-ZB397)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human IL-6
Target	IL6
Immunogen	Recombinant human IL-6
Isotype	IgG2b, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-F5
Conjugate	Unconjugated
Applications	ELISA(Det), ELISPOT Recommended Working Concentration: ELISA / ELISpot: 0.05-0.5 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB394 - CABT-ZB397
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS
Preservative	None

Storage	Stable at +2-8°C for 12 months
----------------	--------------------------------

Ship	Wet ice
-------------	---------

GENE INFORMATION

Synonyms	IL6; interleukin 6; HGF; HSF; BSF2; IL-6; IFNB2; interleukin-6; CDF; BSF-2; IFN-beta-2; interferon beta-2; interleukin BSF-2; interferon; beta 2; hybridoma growth factor; CTL differentiation factor; B-cell stimulatory factor 2; B-cell differentiation factor; PF-04236921
-----------------	--

Entrez Gene ID	3569
-----------------------	----------------------

UniProt ID	P05231
-------------------	------------------------