



# Mouse Anti-Human IL-7 monoclonal antibody, clone C-O29 (CABT-ZB399)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Recognizes both natural and recombinant human IL-7
<b>Target</b>	IL7
<b>Immunogen</b>	Recombinant human IL-7
<b>Isotype</b>	IgG1, κ
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	C-O29
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA(Cap) Recommended Working Concentration: ELISA: 0.5-2.5 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB399 - CABT-ZB401
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1.0 mg
<b>Buffer</b>	PBS
<b>Preservative</b>	None

<b>Storage</b>	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
----------------	---

<b>Ship</b>	Wet ice
-------------	---------

## GENE INFORMATION

<b>Synonyms</b>	IL7; interleukin 7; IL-7; interleukin-7
-----------------	---

<b>Entrez Gene ID</b>	<a href="#">3574</a>
-----------------------	----------------------

<b>UniProt ID</b>	<a href="#">P13232</a>
-------------------	------------------------