



# Mouse Anti-Human IP-10 monoclonal antibody, clone C-D61 (CABT-ZB403)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Recognizes both natural and recombinant human IP-10
<b>Target</b>	CXCL10
<b>Immunogen</b>	Recombinant human IP-10
<b>Isotype</b>	IgG2a, κ
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	C-D61
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	<p>ELISA(Cap), FC</p> <p>Recommended Working Concentration:</p> <p>ELISA: 0.1-1 µg/ml</p> <p>FC: Use 10 µl to label 10<sup>6</sup> cells or 100 µl of whole blood</p> <p>We recommend the following for sandwich ELISA (Capture - Detection):</p> <p>CABT-ZB403 - CABT-ZB407</p> <p>CABT-ZB405 - CABT-ZB407</p>
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1.0 mg
<b>Buffer</b>	PBS

<b>Preservative</b>	None
<b>Storage</b>	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
<b>Ship</b>	Wet ice

## GENE INFORMATION

<b>Synonyms</b>	CXCL10; chemokine (C-X-C motif) ligand 10; C7; IFI10; INP10; IP-10; crg-2; mob-1; SCYB10; gIP-10; C-X-C motif chemokine 10; gamma IP10; gamma-IP10; small-inducible cytokine B10; interferon-inducible cytokine IP-10; 10 kDa interferon gamma-induced protein; protein 10 from interferon (gamma)-induced cell line; small inducible cytokine subfamily B (Cys-X-Cys); member 10; NI-0801; HuAIP12
<b>Entrez Gene ID</b>	<a href="#">3627</a>
<b>UniProt ID</b>	<a href="#">P02778</a>