



Mouse Anti-Human Perforin monoclonal antibody, clone C-F59 (CABT-ZB413)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes human Perforin, a 70 kDa protein
Target	PRF1
Immunogen	Natural human perforin
Isotype	IgG2b, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-F59
Conjugate	Unconjugated
Applications	ELISA(Cap), ELISPOT Recommended Working Concentration: ELISA: 0.5-2.5 µg/ml ELISpot: 5-10 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB413 - CABT-ZB415; CABT-ZB413 - CABT-ZB417; CABT-ZB413 - CABT-ZB418; CABT-ZB413 - CABT-ZB419
Format	Purified, Liquid
Concentration	Lot specific
Size	1.0 mg
Buffer	PBS

Preservative	None
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Ship	Wet ice

GENE INFORMATION

Synonyms	PRF1; perforin 1 (pore forming protein); P1; PFP; FLH2; PFN1; HPLH2; perforin-1; cytolysin; lymphocyte pore forming protein; lymphocyte pore-forming protein
Entrez Gene ID	5551
UniProt ID	P14222