



Mouse Anti-Human TNF- α monoclonal antibody, clone C-G8 (CABT-ZB421)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human TNF- α
Target	TNF
Immunogen	Recombinant human TNF- α
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-G8
Conjugate	Unconjugated
Applications	ELISA(Cap), ELISPOT Recommended Working Concentration: ELISA: 1-10 μ g/ml ELISpot: 5-10 μ g/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB421 - CABT-ZB423
Format	Purified, Liquid
Concentration	Lot specific
Size	1.0 mg
Buffer	PBS

Preservative	None
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Ship	Wet ice

GENE INFORMATION

Synonyms	TNF; DIF; TNFA; TNFSF2; TNF-alpha; TNLG1F; Tumor necrosis factor
Entrez Gene ID	7124
UniProt ID	P01375