



Mouse Anti-Human TNF- α monoclonal antibody, clone C-D8 (CABT-ZB423)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human TNF- α
Target	TNF
Immunogen	Recombinant human TNF- α
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-D8
Conjugate	Unconjugated
Applications	ELISA(Det), ELISPOT, IHC, FA Recommended Working Concentration: ELISA / ELISpot: 0.05-0.5 μ g/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB421 - CABT-ZB423
Format	Purified, Liquid
Concentration	Lot specific
Size	100 μ g
Buffer	PBS
Preservative	None

Storage	Stable at +2-8°C for 12 months
----------------	--------------------------------

Ship	Wet ice
-------------	---------

GENE INFORMATION

Synonyms	TNF; DIF; TNFA; TNFSF2; TNF-alpha; TNLG1F; Tumor necrosis factor
-----------------	--

Entrez Gene ID	7124
-----------------------	----------------------

UniProt ID	P01375
-------------------	------------------------