



Rabbit Anti-FAM polyclonal antibody (CABT-ZB427)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	FAM
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µL
Buffer	PBS
Preservative	None
Storage	Store at 4°C.
Ship	Wet ice

BACKGROUND

Introduction	6-Fam (6-carboxyfluorescein) is a single isomer of the dye Fluorescein. It is usually coupled to DNA and RNA oligonucleotides via an amide bond and has good water solubility. 6-Fam has
--------------	--

the same spectral characteristics as FITC (Fluorescein-5,6-isothiocyanate), which is composed of reactive species of fluorescein with isothiocyanate reactive groups.

Keywords	6-Carboxyfluorescein; 6-FAM; FAM
-----------------	----------------------------------

GENE INFORMATION

Synonyms	6-Carboxyfluorescein; 6-FAM; FAM
-----------------	----------------------------------
