



Mouse Anti-Human MMP-3 monoclonal antibody, clone NN13 (CABT-ZB521)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human MMP-3
Target	MMP3
Immunogen	Recombinant Human MMP-3 Protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	NN13
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(cap) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB521 - CABT-ZB891 This antibody will detect MMP-3 in antibody pair set. [ABPR-ZB097]
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human MMP-3. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Format	Purified, Liquid
Concentration	Lot specific

Size	50 µL, 100 µL, 200 µL, 1 mL
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction Matrix metalloproteinase 3 (abbreviated as MMP3) is also known as stromelysin 1 and progelatinase. MMP3 is a member of the matrix metalloproteinase (MMP) family whose members are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, tissue remodeling, and disease processes including arthritis and metastasis. As a secreted zinc-dependent endopeptidase, MMP3 exerts its functions mainly in the extracellular matrix. This protein is activated by two major endogenous inhibitors: alpha2-macroglobulin and tissue inhibitors of metalloproteinases (TIMPs). MMP3 plays a central role in degrading collagen types II, III, IV, IX, and X, proteoglycans, fibronectin, laminin, and elastin. Also, MMP3 can activate other MMPs such as MMP1, MMP7, and MMP9, rendering MMP3 crucial in connective tissue remodeling. Dysregulation of MMPs has been implicated in many diseases including arthritis, chronic ulcers, encephalomyelitis, and cancer. Synthetic or natural inhibitors of MMPs result in inhibition of metastasis, while up-regulation of MMPs led to enhanced cancer cell invasion.

Keywords MMP3; matrix metalloproteinase 3 (stromelysin 1, progelatinase); SL-1; STMY

GENE INFORMATION

Synonyms MMP3; matrix metalloproteinase 3 (stromelysin 1, progelatinase); SL-1; STMY; STR1; CHDS6; MMP-3; STMY1; stromelysin-1; transin-1; proteoglycanase; matrix metalloproteinase-3; matrix metalloproteinase 3 (stromelysin 1, progelatinase)

Entrez Gene ID [4314](#)

UniProt ID [P08254](#)