



Mouse Anti-Human PBK/TOPK monoclonal antibody, clone NN21 (CABT-ZB526)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human PBK/TOPK
Target	PBK
Immunogen	Recombinant Human PBK/TOPK Protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	NN21
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(cap) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB526 - CABT-ZB893 This antibody will detect PBK/TOPK in antibody pair set. [ABPR-ZB102]
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human PBK / TOPK. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Format	Purified, Liquid
Concentration	Lot specific

Size	50 µL, 100 µL, 200 µL, 1 mL
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	PDZ binding kinase (PBK), also known as TOPK (T-LAK cell-originated protein kinase), is a serine/threonine kinase related to the dual specific mitogen-activated protein kinase kinase (MAPKK) family, and has all the characteristic protein kinase subdomains and a C-terminal PDZ-binding T/SXV motif. PBK is expressed in the testis restrictedly expressed in outer cell layer of seminiferous tubules, as well as placenta. PBK may be enrolled in the activation of lymphoid cells and support testicular functions, with a suggested role in the process of spermatogenesis. This mitotic kinase phosphorylates MAP kinase p38 and seems to be active in mitosis. When phosphorylated, PBK forms a protein-protein interaction with tumor suppressor p53 (TP53), leading to TP53 destabilization and attenuation of G2/M checkpoint during doxorubicin-induced DNA damage. The expression level of PBK is thus upregulated in a variety of neoplasms including hematological malignancies.
Keywords	PBK; PDZ binding kinase; lymphokine-activated killer T-cell-originated protein kinase; cancer/testis antigen 84

GENE INFORMATION

Synonyms	PBK; PDZ binding kinase; lymphokine-activated killer T-cell-originated protein kinase; cancer/testis antigen 84; CT84; FLJ14385; Nori 3; SPK; T LAK cell originated protein kinase; TOPK
Entrez Gene ID	55872
UniProt ID	Q96KB5