



Rabbit Anti-Mouse AXL monoclonal antibody, clone S210 (CABT-ZB554)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Mouse AXL
Target	AXL
Immunogen	Recombinant Mouse AXL Protein
Isotype	IgG1
Source/Host	Rabbit
Species Reactivity	Mouse
Clone	S210
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(cap) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB554 - CABT-ZB916 This antibody will detect AXL in antibody pair set. [ABPR-ZB130]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse AXL.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL

Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Axl receptor tyrosine kinase, together with Tyro3 and Mer, constitute the TAM family of receptor tyrosine kinases. In the nervous system, Axl and its ligand Growth-arrest-specific protein 6 (Gas6) are expressed on multiple cell types. Axl functions in dampening the immune response, regulating cytokine secretion, clearing apoptotic cells and debris, and maintaining cell survival. Axl is upregulated in various disease states, such as in the cuprizone toxicity-induced model of demyelination and in multiple sclerosis (MS) lesions, suggesting that it plays a role in disease pathogenesis. Axl expression correlates with poor prognosis in several cancers. Axl mediates multiple oncogenic phenotypes and activation of these RTKs constitutes a mechanism of chemoresistance in a variety of solid tumors. Axl contributes to cell survival, migration, invasion, metastasis and chemosensitivity justify further investigation of Axl as novel therapeutic targets in cancer. The receptor tyrosine kinase AXL is thought to play a role in metastasis. The soluble AXL receptor as a therapeutic candidate agent for treatment of metastatic ovarian cancer. GAS6/AXL targeting as an effective strategy for inhibition of metastatic tumor progression in vivo.
---------------------	--

Keywords	AXL; AXL receptor tyrosine kinase; UFO; JTK11
-----------------	---

GENE INFORMATION

Synonyms	AXL; AXL receptor tyrosine kinase; UFO; JTK11; tyrosine-protein kinase receptor UFO; AXL oncogene; AXL transforming sequence/gene
-----------------	---

Entrez Gene ID	26362
-----------------------	-----------------------

UniProt ID	Q6PE80
-------------------	------------------------