



Rabbit Anti-Mouse CLEC4N monoclonal antibody, clone S114 (CABT-ZB563)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Mouse CLEC4N
Target	CLEC4N
Immunogen	Recombinant Mouse CLEC4N/CLEC6A/Dectin-2 Protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Clone	S114
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(cap) This antibody will detect CLEC4N in antibody pair set. [ABPR-ZB140]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CLEC4N / CLEC6A / Dectin-2.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL
Buffer	PBS

Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Keywords	CLEC4N; C-type lectin domain family 4, member n; Nkcl; Clec6a
-----------------	---

GENE INFORMATION

Synonyms	CLEC4N; C-type lectin domain family 4, member n; Nkcl; Clec6a; Clecsf10; C-type lectin domain family 6 member A; C-type (calcium dependent, carbohydrate recognition domain) lectin, superfamily member 10; C-type lectin superfamily member 10; DC-associated C-type lectin 2; dectin-2
Entrez Gene ID	56620
UniProt ID	Q9JKF4