



Rabbit Anti-Human IL-18R α monoclonal antibody, clone S122 (CABT-ZB625)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human IL-18R α It has no cross-reactivity in ELISA with Mouse IL18R1/IL1RRP/CD218a, Mouse SIGIRR, Human IL1R1/CD121a, Human IL1R2/IL1RB/CD121b, Human IL1R3/IL1RAP, Human IL1R8/IL1RAPL1, Human IL1RL1/IL1R4/ST2, Human IL18RAP/IL1R7, Human IL1R9/IL1RAPL2.
Target	IL18R1
Immunogen	Recombinant Human IL18R1 protein
Isotype	IgG1
Source/Host	Rabbit
Species Reactivity	Human
Clone	S122
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, ELISA(cap) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB625 - CABT-ZB972 This antibody will detect IL-18R α in antibody pair set. [ABPR-ZB204]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human IL18R1 / CD218a.
Format	Purified, Liquid

Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Interleukin-18 receptor 1 (IL18R1) also known as CD218 antigen-like family member A, CDw218a, IL1 receptor-related protein, and CD218a, is an interleukin receptor of the immunoglobulin superfamily. IL18R1 is found expressed in the lung, leukocytes, spleen, liver, thymus, prostate, small intestine, colon, placenta, and heart, and is absent from the brain, skeletal muscle, pancreas, and kidney. A high level of expression is found in Hodgkin disease cell lines. This receptor is specifically bound to interleukin 18 (IL18) and is essential for IL18 mediated signal transduction. IL18R1 contains 3 Ig-like C2-type (immunoglobulin-like) domains and 1 TIR domain. It is a single-pass type I membrane protein. IFN-alpha and IL12 are reported to induce the expression of this receptor in NK and T cells. The increased expression of IL18R1 may contribute pathogenically to disease and is therefore a potential therapeutic target. The absence of a genetic association in the IL18R1 gene itself suggests regulation from other parts of the genome, or as part of the inflammatory cascade in multiple sclerosis without a prime genetic cause.
Keywords	IL18R1; interleukin 18 receptor 1; CD218a; IL18RA

GENE INFORMATION

Synonyms	IL18R1; interleukin 18 receptor 1; CD218a; IL18RA; IL1RRP; CDw218a; IL-1Rrp; interleukin-18 receptor 1; IL-18R1; IL1R-rp; IL-18R-1; IL18Ralpha2; cytokine receptor; IL1 receptor-related protein; CD218 antigen-like family member A
Entrez Gene ID	16182
UniProt ID	Q61098