



# Rabbit Anti-Human HAI-1 monoclonal antibody, clone S127 (CABT-ZB635)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	It reacts with Human HAI-1
<b>Target</b>	SPINT1
<b>Immunogen</b>	Recombinant Human HAI-1/SPINT1 Protein
<b>Isotype</b>	IgG2b
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Clone</b>	S127
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA(cap) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB635 - CABT-ZB981 This antibody will detect HAI-1 in antibody pair set. [ABPR-ZB214]
<b>Preparation</b>	This antibody was obtained from a rabbit immunized with purified, recombinant Human HAI-1/SPINT1.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	50 µL, 100 µL, 1 mL

<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
<b>Ship</b>	Wet ice

## BACKGROUND

**Keywords** SPINT1; serine peptidase inhibitor, Kunitz type 1; HAI; HAI1

## GENE INFORMATION

**Synonyms** SPINT1; serine peptidase inhibitor, Kunitz type 1; HAI; HAI1; MANSC2; kunitz-type protease inhibitor 1; HAI-1; serine protease inhibitor, Kunitz type 1; hepatocyte growth factor activator inhibitor type 1

**Entrez Gene ID** [6692](#)

**UniProt ID** [O43278](#)