



Rabbit Anti-Mouse IL-16 monoclonal antibody, clone S133 (CABT-ZB656)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Mouse IL-16
Target	IL16
Immunogen	Recombinant Mouse IL16 Protein
Isotype	IgG2b
Source/Host	Rabbit
Species Reactivity	Mouse
Clone	S133
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(cap) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB656 - CABT-ZB996 This antibody will detect IL-16 in antibody pair set. [ABPR-ZB235]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse IL16.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL

Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Keywords Interleukin 17; IL-17

GENE INFORMATION

Synonyms Interleukin 17; IL-17

Entrez Gene ID [16171](#)

UniProt ID [Q62386](#)