



Rabbit Anti-Human IL-12 monoclonal antibody, clone S112 (CABT-ZB790)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human IL-12
Target	IL12
Immunogen	Recombinant Human IL12/IL12A & IL12B Protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S112
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(cap) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB790 - CABT-ZB1092 This antibody will detect IL-12 in antibody pair set. [ABPR-ZB370]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human IL12 / IL12A & IL12B.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL

Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Keywords IL-12b; Interleukin 12b; p40; IL-12b

GENE INFORMATION

Synonyms IL-12b; Interleukin 12b; p40; IL-12b; IL12p40; IL-12p40

Entrez Gene ID [3593](#)

UniProt ID [P29460](#)