



Rabbit Anti-Mouse PIGR monoclonal antibody, clone S151 (CABT-ZB804)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Mouse PIGR
Target	PIGR
Immunogen	Recombinant Mouse PIGR protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Clone	S151
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, ELISA(cap) This antibody will detect PIGR in antibody pair set. [ABPR-ZB384]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse PIGR extracellular domain.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL
Buffer	PBS

Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Polymeric immunoglobulin receptor, also known as PIGR, is a member of the immunoglobulin superfamily and a Fc receptor. The ectodomain of this receptor consists of five units with homology to the variable units of immunoglobulins and a transmembrane region, which also has some homology to certain immunoglobulin variable regions. PIGR is expressed on several glandular epithelia including those of liver and breast. The deduced amino-acid sequence has a length of 764 residues and shows an overall similarity of 56% and 64% with the rabbit and rat counterpart. PIGR mediates transcellular transport of polymeric immunoglobulin molecules, and thus facilitates the secretion of IgA and IgM. During this process, a cleavage occurs that separates the extracellular (known as the secretory component) from the transmembrane segment of PIGR.
Keywords	PIGR; polymeric immunoglobulin receptor; poly-Ig receptor; hepatocellular carcinoma associated protein TB6

GENE INFORMATION

Synonyms	PIGR; polymeric immunoglobulin receptor; poly-Ig receptor; hepatocellular carcinoma associated protein TB6; hepatocellular carcinoma-associated protein TB6; FLJ22667; MGC125361; MGC125362
Entrez Gene ID	18703
UniProt ID	O70570