



Rabbit Anti-Rat PCSK9 monoclonal antibody, clone S118 (CABT-ZB864)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Rat PCSK9
Target	PCSK9
Immunogen	Recombinant Rat PCSK9/NARC1 Protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Clone	S118
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB488 - CABT-ZB864 This antibody will detect PCSK9 in antibody pair set. [ABPR-ZB063]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Rat PCSK9 / NARC1.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL

Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Proprotein convertase subtilisin/kexin type 9 (PCSK9), also known as NARC1 (neural apoptosis regulated convertase), which is a newly identified human secretory subtilase belonging to the proteinase K subfamily of the secretory subtilase family. PCSK9 protein is an enzyme which in humans is encoded by the PCSK9 gene with orthologs found across many species. It is expressed in neuroepithelioma, colon carcinoma, hepatic and pancreatic cell lines, and in Schwann cells. PCSK9 protein is highly expressed in the liver and regulates low density lipoprotein receptor (LDLR) protein levels. Inhibition of PCSK9 protein function is currently being explored as a means of lowering cholesterol levels. Thereby, PCSK9 protein is regarded as a new strategy to treat hypercholesterolemia. PCSK9 protein contributes to cholesterol homeostasis and may have a role in the differentiation of cortical neurons.
---------------------	---

Keywords	PCSK9; proprotein convertase subtilisin/kexin type 9; FH3; PC9
-----------------	--

GENE INFORMATION

Synonyms	PCSK9; proprotein convertase subtilisin/kexin type 9; FH3; PC9; NARC1; LDLCQ1; NARC-1; HCHOLA3; subtilisin/kexin-like protease PC9; neural apoptosis regulated convertase 1; convertase subtilisin/kexin type 9 preproprotein; ALD306
-----------------	---

Entrez Gene ID	298296
-----------------------	------------------------

UniProt ID	P59996
-------------------	------------------------