



Rabbit Anti-Human CD32B/Fcgr2b monoclonal antibody, clone S223 (CABT-ZB866)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	<p>It reacts with Human CD32B/Fcgr2b</p> <p>It has cross-reactivity in ELISA with Human CD64/FCGR1A.</p> <p>It has no cross-reactivity in ELISA with Mouse CD32b/FCGR2B, Human CD32a/FCGR2A, Human CD16b/FCGR3B, Human CD16a/FCGR3A.</p>
Target	FCGR2B
Immunogen	Recombinant Human CD32b protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S223
Purification	Protein A purified
Conjugate	Unconjugated
Applications	<p>ELISA, ELISA(det)</p> <p>We recommend the following for sandwich ELISA (Capture - Detection):</p> <p>CABT-ZB492 - CABT-ZB866</p> <p>This antibody will detect CD32B/Fcgr2b in antibody pair set. [ABPR-ZB067]</p>
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CD32b extracellular domain.
Format	Purified, Liquid

Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Keywords FCGR2B; Fc fragment of IgG, low affinity IIb, receptor (CD32); CD32; FCG2

GENE INFORMATION

Synonyms	FCGR2B; Fc fragment of IgG, low affinity IIb, receptor (CD32); CD32; FCG2; CD32B; FCGR2; IGFR2; low affinity immunoglobulin gamma Fc region receptor II-b; CDw32; fcRII-b; Fc gamma RIIB; fc-gamma-RIIb; fc-gamma RII-b; igG Fc receptor II-b; Fc fragment of IgG, low affinity II, receptor for (CD32); Fc fragment of IgG, low affinity IIb, receptor for (CD32)
Entrez Gene ID	2213
UniProt ID	P31994