



Rabbit Anti-Human TNF beta monoclonal antibody, clone S185 (CABT-ZB885)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human TNF beta It has no cross-reactivity in ELISA with Human TNFa, Mouse TNFa, Human TNFRII.
Target	TNF
Immunogen	Recombinant Human TNF-beta/TNFSF1/Lymphotoxin alpha protein
Isotype	IgG2a
Source/Host	Rabbit
Species Reactivity	Human
Clone	S185
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB515 - CABT-ZB885 This antibody will detect TNF beta in antibody pair set. [ABPR-ZB091]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human TNF-beta / TNFSF1 / Lymphotoxin alpha .
Format	Purified, Liquid
Concentration	Lot specific

Size	50 µL, 100 µL, 1 mL
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	<p>Lymphotoxin-alpha, also known as LT-alpha, TNF-beta, Tumor necrosis factor ligand superfamily member 1, LTA TNFSF1, and TNFB, is a secreted protein that belongs to the tumor necrosis factor family. TNF-beta/TNFSF1/Lymphotoxin alpha is highly inducible, secreted, and exists as a homotrimeric molecule. It is a cytokine that in its homotrimeric form binds to TNFRSF1A/TNFR1, TNFRSF1B/TNFR1, and TNFRSF14/HVEM. In its heterotrimeric form with LTB, TNF-beta/TNFSF1/Lymphotoxin alpha binds to TNFRSF3/LTBR. Lymphotoxin is produced by lymphocytes and cytotoxic for a wide range of tumor cells. TNF-beta/TNFSF1/Lymphotoxin alpha forms heterotrimers with lymphotoxin-beta which anchors lymphotoxin-alpha to the cell surface. It mediates a large variety of inflammatory, immunostimulatory, and antiviral responses. TNF-beta/TNFSF1/Lymphotoxin alpha is also involved in the formation of secondary lymphoid organs during development and plays a role in apoptosis. Genetic variations in TNF-beta/TNFSF1/Lymphotoxin alpha are a cause of susceptibility psoriatic arthritis which is an inflammatory, seronegative arthritis associated with psoriasis. It is a heterogeneous disorder ranging from a mild, non-destructive disease to a severe, progressive, erosive arthropathy.</p>
---------------------	---

Keywords	TNF; DIF; TNFA; TNFSF2
-----------------	------------------------

GENE INFORMATION

Synonyms	TNF; DIF; TNFA; TNFSF2; TNF-alpha; TNLG1F; Tumor necrosis factor
Entrez Gene ID	7124
UniProt ID	Q9UNG2