



Rabbit Anti-Human SAA4 monoclonal antibody, clone S125 (CABT-ZB898)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human SAA4
Target	SAA4
Immunogen	Recombinant Human SAA4 Protein
Isotype	IgG1
Source/Host	Rabbit
Species Reactivity	Human
Clone	S125
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB534 - CABT-ZB898 This antibody will detect SAA4 in antibody pair set. [ABPR-ZB110]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human SAA4.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL

Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	SAA4 is a member of the SAA family. SAA proteins are family of apolipoproteins of high density lipoprotein (HDL). They can be separated into two distinct groups. First group (SAA1, SAA2, and SAA3) consists of acute phase reactant whose expression level increase in the blood in a response to trauma, infection, inflammation, and neoplasia. These acute phase SAAs associates with HDL during inflammation and remodel the HDL particle by displacing Apo-A1. The second distinct group consists of SAA4 and SAA5 which exist as the minor apolipoproteins on HDL, but this group of SAA constitutes more than 90% of all the SAA during homeostasis, and it is thought to play a role in the normal functioning of the HDL particle. SAA4 is a constitutively expressed protein which expressed only in humans and mice. It is connected almost completely with lipoproteins of the high density range. The physiological function of SAA4 is unknown, and its serum concentration has no association with those of other major apolipoproteins.
---------------------	---

Keywords	SAA4; serum amyloid A4, constitutive; serum amyloid A-4 protein; C SAA
-----------------	--

GENE INFORMATION

Synonyms	SAA4; serum amyloid A4, constitutive; serum amyloid A-4 protein; C SAA; CSAA; constitutively expressed serum amyloid A protein; C-SAA
-----------------	---

Entrez Gene ID	6291
-----------------------	----------------------

UniProt ID	P35542
-------------------	------------------------