



Mouse Anti-Rhesus Monkey CD32A monoclonal antibody, clone NN13 (CABT- ZB927)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Rhesus Monkey CD32A
Target	FCGR2A
Immunogen	Recombinant Rhesus Monkey CD32a/FCGR2A/Fc gamma RIIA Protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Rhesus Monkey
Clone	NN13
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB572 - CABT-ZB927 This antibody will detect CD32A in antibody pair set. [ABPR-ZB149]
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Rhesus CD32a / FCGR2A / Fc gamma RIIA. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Format	Purified, Liquid

Concentration	Lot specific
Size	50 µL, 100 µL, 200 µL, 1 mL
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Keywords FCGR2A; Fc fragment of IgG, low affinity IIa, receptor (CD32); CD32; FCG2

GENE INFORMATION

Synonyms FCGR2A; Fc fragment of IgG, low affinity IIa, receptor (CD32); CD32; FCG2; FcGR; CD32A; CDw32; FCGR2; IGFR2; FCGR2A1; low affinity immunoglobulin gamma Fc region receptor II-a; fcRII-a; fc-gamma-RIIa; fc-gamma RII-a; igG Fc receptor II-a; Immunoglobulin G Fc receptor II