



Rabbit Anti-Rat Beta-2 microglobulin monoclonal antibody, clone S146 (CABT- ZB928)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Rat Beta-2 microglobulin
Target	B2M
Immunogen	Recombinant Rat B2M/beta-2 microglobulin Protein
Isotype	IgG1
Source/Host	Rabbit
Species Reactivity	Rat
Clone	S146
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, ELISA(det), IHC-P We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB573 - CABT-ZB928 This antibody will detect Beta-2 microglobulin in antibody pair set. [ABPR-ZB150]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Rat B2M / beta-2 microglobulin.
Format	Purified, Liquid
Concentration	Lot specific

Size	50 µL, 100 µL, 1 mL
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	B2M, also known as β 2-Microglobulin or CDABP0092, is a component of MHC class I molecules found expression in all nucleated cells (excludes red blood cells). The major function of MHC class I molecules is to display fragments of proteins from within the cell to T-cells and cells containing foreign proteins will be attacked. B2M (β 2-Microglobulin) is a low molecular weight protein. It was demonstrated that B2M (β 2-Microglobulin) was localized in the membranes of nucleated cells and was found to be associated with HL-A antigens. B2M (β 2-Microglobulin) is present in free form in various body fluids and as a subunit of histocompatibility antigens on cell surfaces lateral to the α 3 chain. Unlike α 3, β 2 has no transmembrane region. Directly above β 2 lies the α 1 chain, which itself is lateral to the α 2. In the absence of B2M (β 2 microglobulin), very limited amounts of MHC class I (classical and non-classical) molecules can be detected on the surface. In the absence of MHC class I, CD8 T cells, a subset of T cells involved in the development of acquired immunity cannot develop. Low levels of B2M (β 2 microglobulin) can indicate non-progression of HIV.
---------------------	---

Keywords	B2M; beta-2-microglobulin; Beta-2-microglobulin; beta chain of MHC class I molecules
-----------------	--

GENE INFORMATION

Synonyms	B2M; beta-2-microglobulin; Beta-2-microglobulin; beta chain of MHC class I molecules; beta-2-microglobulin; β 2-MG
Entrez Gene ID	24223
UniProt ID	P07151