



# Rabbit Anti-Human IL-16 monoclonal antibody, clone S339 (CABT-ZB934)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	It reacts with Human IL-16
<b>Target</b>	IL16
<b>Immunogen</b>	Recombinant Human IL16 Protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Clone</b>	S339
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA(det) This antibody will detect in antibody pair set. [ABPR-ZB158]
<b>Preparation</b>	This antibody was obtained from a rabbit immunized with purified, recombinant Human IL16.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	50 µL, 100 µL, 1 mL
<b>Buffer</b>	PBS

<b>Preservative</b>	None
<b>Storage</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
<b>Ship</b>	Wet ice

## BACKGROUND

<b>Keywords</b>	IL16; interleukin 16; LCF; NIL16
-----------------	----------------------------------

## GENE INFORMATION

<b>Synonyms</b>	IL16; interleukin 16; LCF; NIL16; PRIL16; prIL-16; pro-interleukin-16; prointerleukin 16; neuronal interleukin 16; lymphocyte chemoattractant factor
<b>Entrez Gene ID</b>	<a href="#">16170</a>
<b>UniProt ID</b>	<a href="#">O54824</a>