



Mouse Anti-Human IGFBP7 monoclonal antibody, clone NN12 (CABT-ZB935)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human IGFBP7 It has no cross-reactivity in ELISA with Human IGFBP3, Human IGFBP4, Human IGFBP6, Human cell lysate (293 cell line).
Target	IGFBP7
Immunogen	Recombinant Human IGFBP7 protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	NN12
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, ELISA(det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB582 - CABT-ZB935 This antibody will detect IGFBP7 in antibody pair set. [ABPR-ZB159]
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human IGFBP7. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Format	Purified, Liquid

Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Insulin-like growth factor-binding protein 7 (IGFBP7) is a member of the IGFBP family. It has been identified in colorectal adenocarcinoma (CRC) cell lines. The Insulin-like growth factor-binding protein also known as IGFBP serves as a carrier protein for Insulin-like growth factor 1. IGFBPs are distinct but are sharing regions with strong homology. All members of the IGFBP family bind IGF-I and IGF-II with about equal affinity. Insulin-like growth factor (IGF) binding proteins (IGFBPs) have been shown to either inhibit or enhance the action of IGF or act in an IGF-independent manner in the prostate. IGFBP7 could inhibit cell growth, decrease soft agar colony formation activity, and induce apoptosis in RKO and SW620 cells. There is mounting evidence that the structure of the IGFBP proteins plays a key role in the regulation of IGF bioavailability, by modulating its molecular size, capillary membrane permeability, target tissue specificity, cell membrane adherence, and IGF affinity.
Keywords	IL10; interleukin 10; IL10X; interleukin-10

GENE INFORMATION

Synonyms	IL10; interleukin 10; IL10X; interleukin-10; CSIF; IL-10; cytokine synthesis inhibitory factor
Entrez Gene ID	3586
UniProt ID	P22301