



Rabbit Anti-Human ECM1 monoclonal antibody, clone S141 (CABT-ZB955)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human ECM1
Target	ECM1
Immunogen	Recombinant Human ECM1 Protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S141
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB604 - CABT-ZB955 This antibody will detect ECM1 in antibody pair set. [ABPR-ZB182]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human ECM1 .
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL

Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Extracellular matrix protein 1 (ECM1) is a secreted glycoprotein and playing a pivotal role in endochondral bone formation, angiogenesis, and tumour biology. Three splice variants have been identified: ECM1a (540 aa) is most widely expressed, with highest expression in the placenta and heart; ECM1b (415 aa) is differentiation-dependent expressed and found only in tonsil and associated with suprabasal keratinocytes; ECM1c (559 aa) accounts for approximately 15% of skin ECM1. Although ECM1 is not tumor specific, is significantly elevated in many malignant epithelial tumors and is suggested as a possible trigger for angiogenesis, tumor progression and malignancies. It also has been shown to regulate endochondral bone formation, skeletal development and tissue remodeling.
Keywords	ECM1; extracellular matrix protein 1; URBWD; secretory component p85

GENE INFORMATION

Synonyms	ECM1; extracellular matrix protein 1; URBWD; secretory component p85
Entrez Gene ID	1893
UniProt ID	Q16610