



# Rabbit Anti-Human TFPI monoclonal antibody, clone S512 (CABT-ZB963)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	<p>It reacts with Human TFPI</p> <p>It has cross-reactivity in ELISA with Mouse TFPI.</p> <p>It has no cross-reactivity in ELISA with Human TFPI2, Human Vitronectin.</p>
<b>Target</b>	TFPI
<b>Immunogen</b>	Recombinant Human TFPI protein
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Clone</b>	S512
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	<p>ELISA(det), IHC-P, ICC/IF</p> <p>We recommend the following for sandwich ELISA (Capture - Detection):</p> <p>CABT-ZB613 - CABT-ZB963</p> <p>This antibody will detect TFPI in antibody pair set. [ABPR-ZB191]</p>
<b>Preparation</b>	This antibody was obtained from a rabbit immunized with purified, recombinant Human TFPI.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific

<b>Size</b>	50 µL, 100 µL, 1 mL
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
<b>Ship</b>	Wet ice

## BACKGROUND

<b>Introduction</b>	Tissue factor pathway inhibitor (TFPI) is the natural inhibitor of TF coagulant and signaling activities. It is a Kunitz-type serine proteinase inhibitor that down-regulates tissue factor-initiated blood coagulation. With its Kunitz domains, TFPI exhibits significant homology with human inter-alpha-trypsin inhibitor and bovin basic pancreatic trypsin inhibitor. TFPI is the natural inhibitor of TF coagulant and signaling activities. The importance of TFPI in the regulation of blood coagulation is emphasized by how its activity is modulated in human disease. In a factor (F) Xa-dependent feedback system, TFPI binds directly and inhibits the TF-FVII/FVIIa complex. Normally, TFPI exists in plasma both as a full-length molecule and as variably carboxy-terminal truncated forms. TFPI also circulates in complex with plasma lipoproteins. The levels and the dual inhibitor effect of TFPI on FXa and TF-FVII/FVIIa complex offers insight into the mechanisms of various pathological conditions triggered by TF. TFPI may play an important role in modulating TF-induced thrombogenesis and it may also provide a unique therapeutic approach for prophylaxis and/or treatment of various diseases. In addition, studies have shown that TFPI exhibits antiangiogenic and antimetastatic effects in vitro and in vivo. In animal models of experimental metastasis, both circulating and tumor cell-associated TFPI are shown to significantly reduce tumor cell-induced coagulation activation and lung metastasis.
---------------------	---

<b>Keywords</b>	TFPI; tissue factor pathway inhibitor (lipoprotein-associated coagulation inhibitor); EPI; TFI
-----------------	--

## GENE INFORMATION

<b>Synonyms</b>	TFPI; tissue factor pathway inhibitor (lipoprotein-associated coagulation inhibitor); EPI; TFI; LACI; TFPI1; tissue factor pathway inhibitor; anti-convertin; extrinsic pathway inhibitor
-----------------	---

<b>Entrez Gene ID</b>	<a href="#">7035</a>
-----------------------	----------------------

<b>UniProt ID</b>	<a href="#">P10646</a>
-------------------	------------------------