



Rabbit Anti-Human HAI-1 monoclonal antibody, clone S143 (CABT-ZB981)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human HAI-1
Target	SPINT1
Immunogen	Recombinant Human HAI-1/SPINT1 Protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S143
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB635 - CABT-ZB981 This antibody will detect HAI-1 in antibody pair set. [ABPR-ZB214]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human HAI-1/SPINT1.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL

Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Keywords SPINT1; serine peptidase inhibitor, Kunitz type 1; HAI; HAI1

GENE INFORMATION

Synonyms SPINT1; serine peptidase inhibitor, Kunitz type 1; HAI; HAI1; MANSC2; kunitz-type protease inhibitor 1; HAI-1; serine protease inhibitor, Kunitz type 1; hepatocyte growth factor activator inhibitor type 1

Entrez Gene ID [6692](#)

UniProt ID [O43278](#)