



Rabbit Anti-Rat ACE2 monoclonal antibody, clone S114 (CABT-ZB990)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Rat ACE2 It has cross-reactivity in ELISA with Human ACE2, Mouse ACE2, Rhesus Monkey ACE2.
Target	ACE2
Immunogen	Recombinant Rat ACE2/Angiotensin-Converting Enzyme 2 Protein
Isotype	IgG1
Source/Host	Rabbit
Species Reactivity	Rat
Clone	S114
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB647 - CABT-ZB990 This antibody will detect ACE2 in antibody pair set. [ABPR-ZB226]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Rat ACE2 / Angiotensin-Converting Enzyme 2.
Format	Purified, Liquid
Concentration	Lot specific

Size	50 µL, 100 µL, 1 mL
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction

Angiotensin-converting enzyme 2 (ACE2), a first homolog of ACE, regulates the renin angiotensin system (RAS) by counterbalancing ACE activity. Accumulating evidence in recent years has demonstrated a physiological and pathological role of ACE2 in the cardiovascular, renal and respiratory systems. ACE2 also has an important role in blood pressure control. This enzyme, an homolog of ACE, hydrolyzes angiotensin (Ang) I to produce Ang-(1-9), which is subsequently converted into Ang-(1-7) by a neutral endopeptidase and ACE. ACE2 releases Ang-(1-7) more efficiently than its catalysis of Ang-(1-9) by cleavage of Pro(7)-Phe(8) bound in Ang II. Thus, the major biologically active product of ACE2 is Ang-(1-7), which is considered to be a beneficial peptide of the RAS cascade in the cardiovascular system. A physiological role for ACE2 has been implicated in hypertension, cardiac function, heart function and diabetes, and as a receptor of the severe acute respiratory syndrome coronavirus. In the acute respiratory distress syndrome (ARDS), ACE, AngII, and AT1R promote the disease pathogenesis, whereas ACE2 and the AT2R protect from ARDS. Importantly, ACE2 has been identified as a key SARS-coronavirus receptor and plays a protective role in severe acute respiratory syndrome (SARS) pathogenesis. Furthermore, the recent explosion of research into the ACE2 homolog, collectrin, has revealed a new physiological function of ACE2 as an amino acid transporter, which explains the pathogenic role of gene mutations in Hartnup disorder. This review summarizes and discusses the recently unveiled roles for ACE2 in disease pathogenesis.

Keywords	ACE2; angiotensin I converting enzyme 2; ACEH; angiotensin-converting enzyme 2
-----------------	--

GENE INFORMATION

Synonyms

ACE2; angiotensin I converting enzyme 2; ACEH; angiotensin-converting enzyme 2; peptidyl-dipeptidase A; metalloprotease MPROT15; ACE-related carboxypeptidase; angiotensin-converting enzyme homolog; angiotensin I converting enzyme (peptidyl-dipeptidase A) 2

Entrez Gene ID	302668
-----------------------	------------------------

UniProt ID

[Q5EGZ1](#)
