



Anti-ERG monoclonal antibody, clone 6G23 (CABT-ZC1132)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Monoclonal Antibody to Human ERG molecule
Antigen Description	The human paired box protein Pax-5 (PAX5) is a member of the paired box domain family of transcription factors that play a role in B-cell differentiation, neural development, and spermatogenesis. As a result, its expression is detected in the developing central nervous system (CNS) and testis.
Target	ERG
Immunogen	Full length human recombinant protein of human ERG produced in HEK293T cell.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	6G23
Conjugate	Unconjugated
Applications	ELISA, LMNX
Buffer	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol.
Preservative	0.05% Sodium Azide
Storage	Shipped at 4 °C. Upon delivery store at -20 °C. Dilute in PBS (pH7.3) before use. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

GENE INFORMATION

Entrez Gene ID [2078](#)

UniProt ID [B4DN83](#)
