



Anti-HPGD monoclonal antibody, clone 2D20 (CABT-ZC1169)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Monoclonal Antibody to Human HPGD molecule
Antigen Description	HPGD encodes a member of the short-chain nonmetalloenzyme alcohol dehydrogenase protein family. The encoded enzyme is responsible for the metabolism of prostaglandins, which function in a variety of physiologic and cellular processes such as inflammation. Mutations in this gene result in primary autosomal recessive hypertrophic osteoarthropathy and craniosteoarthropathy. Multiple transcript variants encoding different isoforms have been found for this gene.
Target	HPGD
Immunogen	Full length human recombinant protein of human HPGD produced in HEK293T cell.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2D20
Conjugate	Unconjugated
Applications	ELISA, LMNX
Buffer	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol.
Preservative	0.05% Sodium Azide
Storage	Shipped at 4 °C. Upon delivery store at -20 °C. Dilute in PBS (pH7.3) before use. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

GENE INFORMATION

Entrez Gene ID [3248](#)

UniProt ID [B4DU74](#)
