



Anti-ID3 monoclonal antibody, clone 26F4 (CABT-ZC1175)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Monoclonal Antibody to Human ID3 molecule
Antigen Description	The protein encoded by this gene is a helix-loop-helix (HLH) protein that can form heterodimers with other HLH proteins. However, the encoded protein lacks a basic DNA-binding domain and therefore inhibits the DNA binding of any HLH protein with which it interacts.
Target	ID3
Immunogen	Full length human recombinant protein of human ID3 produced in HEK293T cell.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	26F4
Conjugate	Unconjugated
Applications	ELISA, LMNX
Buffer	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose.
Preservative	0.05% Sodium Azide
Storage	Shipped at 4 °C. Upon delivery store at -20 °C. Dilute in PBS (pH7.3) before use. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

GENE INFORMATION

Entrez Gene ID [3399](#)

UniProt ID [Q02535](#)
