



Anti-STAT4 monoclonal antibody, clone 2C20 (CABT-ZC1303)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Monoclonal Antibody to Human STAT4 molecule
Antigen Description	STAT4 carries out a dual function: signal transduction and activation of transcription. Involved in IL12 signaling.
Target	STAT4
Immunogen	Full length human recombinant protein of human STAT4 produced in HEK293T cell.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2C20
Conjugate	Unconjugated
Applications	ELISA, LMNX
Buffer	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol.
Preservative	0.05% Sodium Azide
Storage	Shipped at 4 °C. Upon delivery store at -20 °C. Dilute in PBS (pH7.3) before use. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

GENE INFORMATION

Entrez Gene ID [6775](#)

UniProt ID [Q14765](#)
