



Magic™ Mouse Anti-M. Catarrhalis MhuA Monoclonal antibody, Clone 3F5-5E5 (CABT- ZS1002)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	MAB 3F5-5E5 is specific for the MhuA surface protein of M. catarrhalis.
Target	M. catarrhalis MhuA
Immunogen	M. catarrhalis strain 25240
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	M. catarrhalis
Clone	3F5-5E5
Conjugate	unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	None
Storage	4°C for short term (1 week), store at -20°C to -80°C for long term(1 year); Avoid repeated freeze-thaw cycles

BACKGROUND

Introduction	Moraxella catarrhalis is a respiratory pathogen responsible for acute bacterial otitis media in children and exacerbation of chronic bronchitis in adults. M. catarrhalis strains are frequently resistant to the bactericidal activity of normal human serum.
Keywords	Moraxella catarrhalis; M. catarrhalis; Moraxella catarrhalis MhuA; M. catarrhalis MhuA