



# Magic™ Mouse Anti-M. Catarrhalis Monoclonal antibody, Clone P107 (CABT-ZS1003)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	M. catarrhalis
Immunogen	M. catarrhalis strain 25240
Isotype	IgG
Source/Host	Mouse
Species Reactivity	M. catarrhalis
Clone	P107
Conjugate	unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	None
Storage	4°C for short term (1 week), store at -20°C to -80°C for long term(1 year); Avoid repeated freeze-thaw cycles

## BACKGROUND

**Introduction** Moraxella catarrhalis is a respiratory pathogen responsible for acute bacterial otitis media in

children and exacerbation of chronic bronchitis in adults. *M. catarrhalis* strains are frequently resistant to the bactericidal activity of normal human serum.

---

<b>Keywords</b>	<i>Moraxella catarrhalis</i> ; <i>M. catarrhalis</i>
-----------------	--

---