



Human ABCF3 blocking peptide (CDBP0274)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-ABCF3 antibody
Antigen Description	ABCF3 (ATP-binding cassette, sub-family F (GCN20), member 3) is a protein-coding gene. Diseases associated with ABCF3 include west nile virus, and major depressive disorder. GO annotations related to this gene include ATPase activity and ATP binding. An important paralog of this gene is ABCF1. Displays an antiviral effect against flaviviruses such as west Nile virus (WNV) in the presence of OAS1B (By similarity).
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ABCF3 ATP-binding cassette, sub-family F (GCN20), member 3 [Homo sapiens]
Official Symbol	ABCF3
Synonyms	ABCF3; ATP-binding cassette, sub-family F (GCN20), member 3; ATP-binding cassette sub-family F member 3; EST201864; FLJ11198;
Entrez Gene ID	55324

mRNA Refseq	NM_018358
Protein Refseq	NP_060828
UniProt ID	Q9NUQ8
Chromosome Location	3q25.1-q25.2
Function	ATP binding; ATPase activity; nucleotide binding;