



Human ABHD6 blocking peptide (CDBP0279)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-ABHD6 antibody
Antigen Description	ABHD6 (abhydrolase domain containing 6) is a protein-coding gene. Diseases associated with ABHD6 include cholesterol ester storage disease, and afferent loop syndrome. GO annotations related to this gene include acylglycerol lipase activity.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ABHD6 abhydrolase domain containing 6 [Homo sapiens (human)]
Official Symbol	ABHD6
Synonyms	ABHD6; abhydrolase domain containing 6; monoacylglycerol lipase ABHD6; lipase protein; 2-arachidonoylglycerol hydrolase; abhydrolase domain-containing protein 6;
Entrez Gene ID	57406
mRNA Refseq	NM_020676.5

Protein Refseq	NP_065727.4
UniProt ID	Q9BV23
Chromosome Location	3p14.3
Pathway	Retrograde endocannabinoid signaling, organism-specific biosystem; Retrograde endocannabinoid signaling, conserved biosystem; triacylglycerol degradation, organism-specific biosystem;
Function	acylglycerol lipase activity;