



Human ADAMTSL5 blocking peptide (CDBP0319)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-ADAMTSL5 antibody
Antigen Description	ADAMTSL5 (ADAMTS-like 5) is a protein-coding gene. Diseases associated with ADAMTSL5 include hypertension. GO annotations related to this gene include heparin binding and metalloendopeptidase activity. An important paralog of this gene is ADAMTSL4.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	ADAMTSL5 ADAMTS-like 5 [Homo sapiens (human)]
Official Symbol	ADAMTSL5

Synonyms	ADAMTSL5; ADAMTS-like 5; THSD6; ADAMTS-like protein 5; ADAMTSL-5; thrombospondin, type I, domain containing 6; thrombospondin type-1 domain-containing protein 6;
Entrez Gene ID	339366
mRNA Refseq	NM_213604.2
Protein Refseq	NP_998769.2
UniProt ID	Q0VD77
Chromosome Location	19p13.3
Function	heparin binding; metalloendopeptidase activity; zinc ion binding;