



Human MLLT10 blocking peptide (CDBP0330)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-MLLT10/AF10 antibody
Antigen Description	This gene encodes a transcription factor and has been identified as a partner gene involved in several chromosomal rearrangements resulting in various leukemias. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2010]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	MLLT10 myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila); translocated to, 10 [Homo sapiens (human)]
Official Symbol	MLLT10
Synonyms	MLLT10; myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila); translocated to, 10; AF10; protein AF-10; type I AF10 protein; type IV AF10 protein; type III AF10 protein; ALL1-fused gene from chromosome 10 protein;

Entrez Gene ID	8028
mRNA Refseq	NM_001195626.1
Protein Refseq	NP_001182555.1
UniProt ID	P55197
Chromosome Location	10p12
Function	DNA binding; protein binding; sequence-specific DNA binding transcription factor activity; zinc ion binding;