



Human APOLD1 blocking peptide (CDBP0450)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-APOLD1/VERGE antibody
Antigen Description	APOLD1 is an endothelial cell early response protein that may play a role in regulation of endothelial cell signaling and vascular function.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	APOLD1 apolipoprotein L domain containing 1 [Homo sapiens]
Official Symbol	APOLD1
Synonyms	APOLD1; apolipoprotein L domain containing 1; apolipoprotein L domain-containing protein 1; DKFZP434F0318; FLJ25138; vascular early response gene protein; VERGE; FLJ95166; DKFZp434F0318;
Entrez Gene ID	81575
mRNA Refseq	NM_001130415

Protein Refseq	NP_001123887
UniProt ID	Q96LR9
Chromosome Location	12p13.2
Function	lipid binding;