



Human apoF blocking peptide (CDBP0452)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Apolipoprotein F antibody
Antigen Description	The product of this gene is one of the minor apolipoproteins found in plasma. This protein forms complexes with lipoproteins and may be involved in transport and/or esterification of cholesterol.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	APOF apolipoprotein F [Homo sapiens]
Official Symbol	apoF
Synonyms	APOF; apolipoprotein F; lipid transfer inhibitor protein; LTIP; Apo-F; MGC22520; DKFZp781G18150;
Entrez Gene ID	319
mRNA Refseq	NM_001638

Protein Refseq	NP_001629
UniProt ID	Q13790
Chromosome Location	12q13
Function	cholesterol binding; lipid transporter activity; receptor binding;