



Mouse Apol6 blocking peptide (CDBP0454)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-apolipoprotein L6 (mouse) antibody
Antigen Description	This gene is a member of the apolipoprotein L gene family. The encoded protein is found in the cytoplasm, where it may affect the movement of lipids or allow the binding of lipids to organelles.
Species	Mouse
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	Apol6 apolipoprotein L 6 [Mus musculus]
Official Symbol	Apol6
Synonyms	APOL6; apolipoprotein L 6; apolipoprotein L6; 2310076O14Rik; Apol6
Entrez Gene ID	71939
mRNA Refseq	NM_001163621

Protein Refseq	NP_001157093
Function	molecular_function;