



Human ARL2 blocking peptide (CDBP0482)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-ARL2 antibody
Antigen Description	This gene encodes a small GTP-binding protein of the RAS superfamily which functions as an ADP-ribosylation factor (ARF). The encoded protein is one of a functionally distinct group of ARF-like genes.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ARL2 ADP-ribosylation factor-like 2 [Homo sapiens]
Official Symbol	ARL2
Synonyms	ARL2; ADP-ribosylation factor-like 2; ADP-ribosylation factor-like protein 2; ARFL2;
Entrez Gene ID	402
mRNA Refseq	NM_001199745

Protein Refseq	NP_001186674
UniProt ID	P36404
Chromosome Location	11q13
Function	GDP binding; GTP binding; GTPase activity; GTPase inhibitor activity; nucleotide binding; protein binding;