



Human ARL6IP1 blocking peptide (CDBP0485)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-ARMER antibody
Antigen Description	ARL6IP1 (ADP-ribosylation factor-like 6 interacting protein 1) is a protein-coding gene. Diseases associated with ARL6IP1 include microphthalmia, and nasopharyngitis. An important paralog of this gene is ENSG00000260342.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	None
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	ARL6IP1 ADP-ribosylation factor-like 6 interacting protein 1 [Homo sapiens (human)]
Official Symbol	ARL6IP1
Synonyms	ARL6IP1; ADP-ribosylation factor-like 6 interacting protein 1; AIP1; ARMER; SPG61; ARL6IP; ADP-ribosylation factor-like protein 6-interacting protein 1; aip-1; ARL-6-interacting protein 1;

apoptotic regulator in the membrane of the endoplasmic reticulum;

Entrez Gene ID	23204
mRNA Refseq	NM_015161.1
Protein Refseq	NP_055976.1
UniProt ID	Q15041
Chromosome Location	16p12-p11.2
Function	protein binding;