



Human ARSH blocking peptide (CDBP0503)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Arylsulfatase H antibody
Antigen Description	Sulfatases, such as ARSH, hydrolyze sulfate esters from sulfated steroids, carbohydrates, proteoglycans, and glycolipids. They are involved in hormone biosynthesis, modulation of cell signaling, and degradation of macromolecules.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ARSH arylsulfatase family, member H [Homo sapiens]
Official Symbol	ARSH
Synonyms	ARSH; arylsulfatase family, member H; arylsulfatase H; ASH; sulfatase;
Entrez Gene ID	347527
mRNA Refseq	NM_001011719

Protein Refseq	NP_001011719
UniProt ID	Q5FYA8
Chromosome Location	Xp22.33
Function	arylsulfatase activity; hydrolase activity; metal ion binding;