



Human ASF1A blocking peptide (CDBP0508)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-ASF1A antibody
Antigen Description	This gene encodes a member of the H3/H4 family of histone chaperone proteins and is similar to the anti-silencing function-1 gene in yeast. The protein is a key component of a histone donor complex that functions in nucleosome assembly. It interacts with histones H3 and H4, and functions together with a chromatin assembly factor during DNA replication and repair.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ASF1A ASF1 anti-silencing function 1 homolog A (S. cerevisiae) [Homo sapiens]
Official Symbol	ASF1A
Synonyms	ASF1A; ASF1 anti-silencing function 1 homolog A (S. cerevisiae); histone chaperone ASF1A; CIA; DKFZP547E2110; hCIA; hAsf1; hAsf1a; CCG1-interacting factor A; anti-silencing function 1A; anti-silencing function protein 1 homolog A; CGI-98; HSPC146;
Entrez Gene ID	25842

mRNA Refseq	NM_014034
Protein Refseq	NP_054753
UniProt ID	Q9Y294
Chromosome Location	6q22.31
Function	chromatin binding; histone binding; protein binding;