



Human ASRGL1 blocking peptide (CDBP0512)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-ASRGL1/ALP antibody
Antigen Description	ASRGL1 (asparaginase like 1) is a protein-coding gene. Diseases associated with ASRGL1 include telogen effluvium, and sagittal sinus thrombosis. GO annotations related to this gene include beta-aspartyl-peptidase activity and N4-(beta-N-acetylglucosaminy)-L-asparaginase activity. An important paralog of this gene is AGA.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ASRGL1[Homo sapiens]
Official Symbol	ASRGL1
Synonyms	ASRGL1; ALP; ALP1; asparaginase like 1 protein; FLJ22316;
Entrez Gene ID	114242
UniProt ID	Q7L266

Chromosome Location 11q12.3
