



Mouse Atad5 blocking peptide (CDBP0513)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Atad5 (mouse) antibody
Antigen Description	ATAD5 (ATPase family, AAA domain containing 5) is a protein-coding gene. Diseases associated with ATAD5 include sotos syndrome, and neurofibromatosis. GO annotations related to this gene include nucleoside-triphosphatase activity and ATP binding.
Species	Mouse
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ATAD5 ATPase family, AAA domain containing 5 [Homo sapiens]
Official Symbol	ATAD5
Synonyms	ATAD5; ATPase family, AAA domain containing 5; C17orf41, chromosome 17 open reading frame 41; ATPase family AAA domain-containing protein 5; ELG1; enhanced level of genomic instability 1 homolog (S. cerevisiae); FLJ12735; FRAG1; chromosome fragility associated gene 1; chromosome fragility-associated gene 1 protein; enhanced level of genomic instability 1 homolog; C17orf41;

Entrez Gene ID	79915
mRNA Refseq	NM_024857
Protein Refseq	NP_079133
UniProt ID	Q96QE3
Chromosome Location	17q11.2
Function	ATP binding; nucleoside-triphosphatase activity; nucleotide binding;