



Human ATG13 blocking peptide (CDBP0519)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-ATG13 antibody
Antigen Description	ATG13 (autophagy related 13) is a protein-coding gene. GO annotations related to this gene include protein kinase binding.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	ATG13 autophagy related 13 [Homo sapiens (human)]
Official Symbol	ATG13
Synonyms	ATG13; autophagy related 13; KIAA0652; PARATARG8; autophagy-related protein 13; ATG13 autophagy related 13 homolog;

Entrez Gene ID	9776
mRNA Refseq	NM_001142673.2
Protein Refseq	NP_001136145.1
UniProt ID	A8K0S6
Chromosome Location	11p11.2
Pathway	Regulation of autophagy, organism-specific biosystem; Regulation of autophagy, conserved biosystem; Senescence and Autophagy, organism-specific biosystem; mTOR signaling pathway, organism-specific biosystem;
Function	protein binding; protein kinase binding;